
## Breathometer Mint The Smart Way to Monitor Oral Health






User Mint Mouthpiece Physical Components Manual Device ..... USB Charging Cord Power **Button** Micro USB\_

**Charging Port** 

## **Smart Components**




Mint test analyzes breath through a highly sophisticated array of sensors, Mint can detect the signs that harmful bacteria leave behind and accurately report through the Mint app. Reports also on breath quality and hydration.



The mint breathometer utilizes an electrochemical printed fuel-cell sensor technology to run its tests.

Identifu Parts



Intake fan (draws in air from mouth to run tests)

## Connectivity Components



Connects to iPhones with iOS 9.1 or newer and Android 4.4.4 or newer via bluetooth LE

